ZIRCONIUM DISCS / BLOCKS

- Isostatically pressed, pre-sintered

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- isostatisch geperst, pre-gesinterd

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- Prensado isostaticamente, pré-sinterizado

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HENRY SCHEIN INC.
(manufacturer)
135 DURYEA ROAD
Melville, NY 11747 USA
EQ REP

HENRY SCHEIN SERVICES GMBH
MONZASTRASSE 2A
D-63225 LANGEN, GERMANY
Rev. 2017/03

0482

Zirlux 16+ pré-pigmentados e brancos

Cronograma de sinterização 2 para Zirlux 16+ pré-pigmentado (usado para pré-pigmentado em branco)				
Procedimento	Temperatura inicial (°C)	Temperatura final (°C)	Tempo de espera (min)	Velocidade de aquecimento (°C/min)
Etapa 1	Temperatura ambiente	1200	0	8
Etapa 2	1200	1200	30	0
Etapa 3	1200	1300	0	2
Etapa 4	1300	1530	0	4
Etapa 5	1530	1530	120	0
Etapa 6	1530	800	0	-8
Etapa 7	800	100	Resfriamento do forno	

Tempo de execução do programa: 495,63 min				
Cronograma de sinterização 3 para Zirlux 16+ pré-pigmentado (para forno de 4 etapas)				
Procedimento	Temperatura inicial (°C)	Temperatura final (°C)	Tempo de espera (min)	Velocidade de aquecimento (°C/min)
Etapa 1	Room Temperature	1200	30	5
Etapa 2	1200	1530	0	2
Etapa 3	1530	1530	120	0
Etapa 4	1530	100	Furnace cooling	
The running time of program : 520min				

Sintering schedule 1 for Zirlux16+ White (used for white blank shaded by coloring liquid)				
Procedure	Start temperature (☒)	Finishing temperature ☒	Holding Time (min)	Heating velocity (☒)/min
Step1	Room Temperature	500	0	2.5
Step2	500	1200	0	8
Step3	1200	1200	30	0
Step4	1200	1300	0	2
Step5	1300	1530	0	4
Step6	1530	1530	120	0
Step7	1530	800	0	-8
Step8	800	100	Furnace cooling	
The running time of program : 626.25min				

Sintering schedule 4 for Zirlux 16+ white with coloring liquids (for 4-step furnace)				
Procedure	Start Temperature (°C)	End Temperature (°C)	Holding Time (min.)	Heat Velocity (°C/min.)
Step 1	Room temperature	500	15	2.5
Step 2	500	1200	30	8
Step 3	1200	1530	120	2
Step 4	1530	100	furnace cooling	

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Zirconium Discs / Blocks

INSTRUCTIONS FOR USE

1. Fix the blank / block according to the CAM-System instructions. Set the enlargement-factor into the unit..
2. Start the milling unit.
3. After milling, extract the blank / block together with the contained milled restoration.
4. Remove the restoration from the blank / block with a fissure bur.
5. Remove any dust from the restoration with compressed air before sintering.
6. Position the restoration on the sinter bed of the furnace tray.
7. Program your furnace according to sintering temperature / instruction shown on the front page.
8. Run the sintering cycle.
9. Wait until cool down is completed before removing the sintered restoration.
10. Inspect the restoration in respect of flaws, wall thickness, seat and marginal fit. If necessary, make small adjustments with a water-cooled diamond bur.
11. Rinse the framework/anatomical crown with water and dry it. It is now ready for veneering. The framework is now ready for veneering /anatomical crown is ready for glazing.!

Safety precautions

- Design the restoration to maximize cross-sections.
- Carefully read the MSDS.

Indications for use

Zirconium discs / blocks are indicated for use as anatomy crown and as substructure for porcelain fused ceramic fixed dental restorations. Limitations are listed in table 1.

Process chain	Material	Single unit crowns		Bridges, number of pontics	
		Anterior	Posterior	Anterior	Posterior
Porous blank/block is machined to enlarged framework/autonomic crown, then sintered and veneered	3Y-TZP	X	X	3	2

*X indicated; digits show maximum number of pontics. Limitation of Bridge, number of Pontics - All Zirlux-discs / blocks are solely for use by or on the order of a dental professional. They are not for use by the general public or over-the-counter.

Contraindications

Contraindications are all applications not indicated in table 1 above.

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